


FossManager™

Secure your instrument configuration and analytic data



ANALYTICS BEYOND MEASURE



Above instruments and equipment are not included.

FossManager™ and **FossManager™ Pro** empowers you to manage your analytical instruments directly from your PC anytime and anywhere.

With the FossManager tool suite, you can configure settings and monitor performance on one or more instruments at a single site or across several locations. Remote management allows you to access your data and back-up settings whenever it suits you and from any location.

If you are concerned about accuracy, compliance and reducing cost of ownership, this tool suite is the perfect solution, allowing you to take full control of your analytical performance.

FossManager and FossManager Pro are supported by Microsoft Azure, safely connecting your instruments via FossConnect.

Your company benefits


- Access instrument settings and consolidated data on a single or multiple instruments across sites, from any location
- Perform remote configuration of your instruments and monitor that they are operated in compliance with standard operating procedures
- Safeguard your data with automatic backups of settings and analytical results

FossManager™

The easiest way to manage your analytical performance!

Feature	FossManager™	FossManager™ Pro
Upload of instrument diagnostics	✓	✓
Automatic backup of sample data	5 years	5 years
Instrument configuration stored in cloud*	✓	✓
Remote network management	✓	✓
Central access to results and data	✓	✓
Access to instrument via FossAssure™ portal	✓	✓
Annual instrument health check		✓
Configuration management by FOSS		✓
FossAPI™ access license	Optional add-on	Included

* FossConnect™ solution supported by Microsoft Azure



FOSS

Tel.: +45 7010 3370

info@foss.dk · www.fossanalytics.com

GB, September 2018