

FossConnect™

Connecting your instruments to a network

	FossConnect™ Basic Powered by MOSAIC	FossConnect™ Manager Powered by MOSAIC	FossConnect™ Manager Plus Powered by MOSAIC
Upload of instrument diagnostics	☑	☑	☑
Automatic backup of sample data	1 YEAR	5 YEARS	5 YEARS
Instrument configuration stored	LOCAL	NETWORK	NETWORK
Remote network management		☑	☑
Central access to results and data		☑	☑
Access to instrument via FossAssure portal		☑	☑
Annual instrument health check			☑
Configuration management service by FOSS			☑
Mosaic API access license*		OPTIONAL ADD-ON	INCLUDED

ANALYTICS BEYOND MEASURE

FossConnect™ can help you reduce your production costs and secure the quality of your products regardless of where they are produced. Besides, FossConnect will also ensure you have peace of mind knowing that your data is backed-up and protected for whenever you need it.

FossConnect™ Basic links your instruments to a secure FOSS server, enabling our technical experts to access vital sample and diagnostics data in order to support you faster and better.

FossConnect™ Manager is a centralised network management solution enabling remote configuration of your entire instrument network while providing a tool for data management and reporting/exporting.

FossConnect™ Manager Plus adds the expert assistance of a FOSS specialist to help you manage your instrument and network. This service also includes a yearly instrument status report and an API license that enables data transfer from Mosaic to LIMS, ERP, and more.


Main Benefits

- Access consolidated data from a centralised location
- Get remote configuration of your instruments in compliance with standard operating procedures
- Secure your data with automatic backups of your configuration and sample data



FossConnect is based on Mosaic, a software platform developed by FOSS, in collaboration with our customers, over the last 10 years.

Connecting your instruments to a network



Instrument availability

All FOSS instruments* can be configured and managed through FossConnect Manager and FossConnect Manager Plus.

FossConnect Basic will only become available for instruments running ISIsScan Nova or Foss Integrator SW.

Currently the following instruments are supported by FossConnect Basic:

- NIRS DS2500 / NIRS DS2500 F
- InfraXact
- Infratec NOVA

*With a few exceptions for older, non-supported instrument models

FOSS
Foss Allé 1
DK-3400 Hilleroed
Denmark

Tel.: +45 7010 3370
Fax: +45 7010 3371

info@foss.dk
www.fossanalytics.com

January 2018