

Recommended PC Requirements for FOSS Software

This document describes the recommended PC requirements for all following FOSS Software for all instruments: Compass, ISIScan Nova, ISIScan, Foss Calibrator, Foss Integrator, Foss Manager and WinISI.

Technical Data

Component	Requirement	Special Instrument Requirements*
Operating Systems	Windows 10 (64 bit) Enterprise, Professional (build 1803 or higher)	
	Windows 10 (64 bit) Enterprise, Windows server 2016 (64 bit)	ProFoss
CPU	Dual-core Intel i5 or equivalent	
CPU speed	2.8 GHz	
Memory (RAM)	8 GB	
	16 GB	Profoss
HDD	SSD disk Minimum 50 GB free HDD space NTFS file system	
	500 GB minimum	ProFoss, CMT instruments
Video	SVGA at 1920x1080, min. 32-bit colour	
NIC	100 Mbit Ethernet IPv4 (IPv6 not supported)	ISIScan Nova
	2 LAN ports	ISIScan Nova CMT instruments (generation 7)
USB	minimum three ports recommended	
COM port	Serial communication port (for use with WineScan Auto Sampler or MilkoScan FT2 balance)	WineScan FT2

*The listed instruments have a higher requirement than the standard recommended requirements.

ANALYTICS BEYOND MEASURE

FOSS Analytical A/S
Nils Foss Allé 1
DK-3400 Hillerød
Denmark
Phone +45 70 10 33 70
Fax +45 70 10 33 71
E-mail info@foss.dk
Web www.fossanalytics.com

FOSS Analytical Co., Ltd.
6 Louyang Road, Building 1
215121, SIP, Suzhou
P.R. China
Phone +86 512 62 92 01 00
Fax +86 512 62 80 56 30
E-mail info@foss.dk
Web www.fossanalytics.com